

Wholly-owned subsidiary Via dei Fornaciai, 9 - Sar 06081 Assisi - Perugia -	y of Sir Holding Srl Viale nta Maria degli Angeli 2009(Italy Tel. +	Nr. IT 03359340548 Europa 74 0 Cusago Milano ·39 02.90394575 (R.A.) ·39 02.90394441	export@sirsafety.com www.sirsafety.com	SATISCATION
Name: Splinter Removal Kit		Code: 81164	First Aid Equipment	
 Description Contains instruments and instructions for removing splinters of various kinds and performing first aid actions. Content: removal pencil mirror magnifying glass eye wash solution 100 ml with an eye glass stainless steel tweezers 11.5 cm packs of sterile gauze 18x40 cm. 				
Category	Colours Blue		Available Sizes Unica	
Quantity/Carton 1	Important Information Read carefully the user sheet included. Make sure the item is in good condition, undamaged before use. Once you become aware of a damage, you should immediately replace it.		Recommended for Suitable for working environments with risk of injured by metal splinters.	
Dexterity	EC Standards		Storage and cleaning Keep in a cool dry place, far from heat sources	6.